



D107A mouse mAb

Catalog No	YP-mAb-09670
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DEFB107A DEFB107 DEFB7; DEFB107B
Immunogen	Synthesized peptide derived from human D107A AA range: 3-53
Specificity	This antibody detects endogenous levels of human D107A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD17 antigen like family member A;CD17 antigen-like family member A;CD17a; CD17a antigen;LAMP 1;LAMP-1;LAMP1;LAMP1_HUMAN;LAMP1;LGP12;lgpA; Lysosomal membrane glycoprotein 12KD;Lysosomal Associated Membrane Protein 1;Lysosome associated membrane glycoprotein 1;Lysosome-associated membrane glycoprotein 1;Lysosome-associated membrane protein 1; OTTHUMP4663;D107A;DEFB107A DEFB107 DEFB7;DEFB107B;Beta-defensin 107 (Beta-defensin 7) (BD-7) (DEFB-7) (Defensin;beta 107); beta 107)
Observed Band	
Calculated Molecular Weight	45kDa,13kDa
Cell Pathway	Secreted.
Tissue Specificity	Specifically expressed in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images