



CINP mouse mAb

Catalog No	YP-mAb-09681
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CINP
Immunogen	Synthesized peptide derived from human CINP AA range: 160-210
Specificity	This antibody detects endogenous levels of human CINP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A disintegrin and metalloproteinase with thrombospondin motifs 2;A disintegrin like and metalloprotease (repolysin type) with thrombospondin type 1 motif 2; ADAM metalloproteinase with thrombospondin type 1 motif 2;ADAM TS 2;ADAM TS2;ADAM-TS 2;ADAM-TS2;ADAMTS 3;ADAMTS-2;ADAMTS2;ATS2_HUMAN; EC 3.4.24.14;EDS VIIB;EDS VIIC;hPCPNI;NPI;PC I NP;PC I-NP;PCINP;PCPNI; pNPI;Procollagen I N proteinase;Procollagen I N-proteinase;Procollagen I/II amino propeptide processing enzyme;Procollagen I/II amino propeptide-processing enzyme;Procollagen N endopeptidase;Procollagen N-endopeptidase; CINP;Cyclin-dependent kinase 2-interacting protein (CDK2-interacting protein)
Observed Band	
Calculated Molecular Weight	135kDa; 135kD
Cell Pathway	Nucleus . Binds to nuclear under G1 conditions, and dissociates from chromatin with the start of DNA replication.
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images