



CPT2 mouse mAb

Catalog No	YP-mAb-09692
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CPT2 CPT1
Immunogen	Synthesized peptide derived from human CPT2 AA range: 459-509
Specificity	This antibody detects endogenous levels of human CPT2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Carnitine O-palmitoyltransferase 2;mitochondrial;carnitine palmitoyltransferase 2; CPT II;CPTASE;EC:2.3.1.21;CPT2;CPT2 CPT1;mitochondrial (EC 2.3.1.21) (Carnitine palmitoyltransferase II) (CPT II); mitochondrial (EC 2.3.1.21) (Carnitine palmitoyltransferase II) (CPT II)
Observed Band	65-70 kDa
Calculated Molecular Weight	74 kDa
Cell Pathway	Mitochondrion inner membrane; Peripheral membrane protein; Matrix side.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images