



SP9 mouse mAb

Catalog No	YP-mAb-09720
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SP9
Immunogen	Synthesized peptide derived from human SP9 AA range: 187-237
Specificity	This antibody detects endogenous levels of human SP9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	bA416N4.2;Bdr;CMS18;dJ168F16.2;FLJ2379;HGNC:11132;MGC15414;MGC139754;Resistance to inhibitors of cholinesterase 4 homolog;RIC 4;RIC4;SEC 9;SEC9;SNAP 25;SNAP;SNAP-25;SNAP-25B;SNAP25;SNP 25;SNP25;SNP25 HUMAN;sp;SUP;Super protein;Synaptosomal associated 25 kDa protein;Synaptosomal associated protein;Synaptosomal associated protein 25;Synaptosomal associated protein 25kDa;Synaptosomal-associated 25 kDa protein;Synaptosomal-associated protein 25;Synaptosomal-associated protein;25-KD;SP9;Transcription factor Sp9
Observed Band	
Calculated Molecular Weight	25kDa; 23kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images