



BDP1 mouse mAb

Catalog No	YP-mAb-09727
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BDP1 KIAA1241 KIAA1689 TFNR
Immunogen	Synthesized peptide derived from human BDP1 AA range: 1732-1782
Specificity	This antibody detects endogenous levels of human BDP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BDP;NEC 1;NEC1;NEC1_HUMAN;Neuroendocrine convertase 1;PC1;PC3;PCSK1;Prohormone convertase 1;Prohormone convertase 1/3;Prohormone convertase 3;Proprotein convertase 1;SPC3;BDP1;BDP1 KIAA1241 KIAA1689 TFNR;Transcription factor TFIIIB component B" homolog (Transcription factor IIB 150) (TFIIIB150) (Transcription factor-like nuclear regulator)
Observed Band	
Calculated Molecular Weight	84kD
Cell Pathway	Nucleus .
Tissue Specificity	Isoform 2 is highly expressed in cerebellum.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images