



NFU1 mouse mAb

Catalog No	YP-mAb-09749
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NFU1 HIRIP5 CGI-33
Immunogen	Synthesized peptide derived from human NFU1 AA range: 203-253
Specificity	This antibody detects endogenous levels of human NFU1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGI 33;HIRA interacting protein 5;HIRIP5;Nfu;NFU1;NifU;NIFUC;NFU1 HIRIP5 CGI-33;NFU1 iron-sulfur cluster scaffold homolog;mitochondrial (HIRA-interacting protein 5); mitochondrial (HIRA-interacting protein 5)
Observed Band	23-28 kDa
Calculated Molecular Weight	254 aa, 28 kDa
Cell Pathway	Mitochondrion. Cytoplasm, cytosol.
Tissue Specificity	Ubiquitous. Expression in adult lung is weak compared to fetal lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images