



RN181 mouse mAb

Catalog No	YP-mAb-09751
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RNF181 HSPC238
Immunogen	Synthesized peptide derived from human RN181 AA range: 7-57
Specificity	This antibody detects endogenous levels of human RN181
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RN181;RING finger protein 181;HSPC238;EC:2.3.2.27;E3 ubiquitin-protein ligase RNF181;RNF181 HSPC238;E3 ubiquitin-protein ligase RNF181 (EC 6.3.2.-) (RING finger protein 181)
Observed Band	18 kDa
Calculated Molecular Weight	153 aa, 18 kDa
Cell Pathway	
Tissue Specificity	Widely expressed, with highest levels in liver and heart and lowest levels in brain and skeletal muscle. Expressed in platelets (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images