



O11H2 mouse mAb

Catalog No	YP-mAb-09752
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	OR11H2 C14orf15 OR11H2P
Immunogen	Synthesized peptide derived from human O11H2 AA range: 269-319
Specificity	This antibody detects endogenous levels of human O11H2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	O11H2;OR11H2 C14orf15 OR11H2P;Olfactory receptor 11H2 (Olfactory receptor OR14-1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images