



LNx1 mouse mAb

Catalog No	YP-mAb-09753
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LNx1 LNx PDZRN2 UNQ574/PRO1136
Immunogen	Synthesized peptide derived from human LNx1 AA range: 527-577
Specificity	This antibody detects endogenous levels of human LNx1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ligand of numb protein X 1;LNx;LNx1;MPDZ;Numb binding protein 1;PDZRN2; LNx1 LNx PDZRN2 UNQ574/PRO1136;E3 ubiquitin-protein ligase LNx (EC 6.3.2.-) (Ligand of Numb-protein X 1) (Numb-binding protein 1) (PDZ domain-containing RING finger protein 2)
Observed Band	70 kDa
Calculated Molecular Weight	728 aa, 81 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in heart, placenta, kidney, pancreas and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images