



ALKB8 mouse mAb

Catalog No	YP-mAb-09754
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ALKBH8 ABH8
Immunogen	Synthesized peptide derived from human ALKB8 AA range: 27-77
Specificity	This antibody detects endogenous levels of human ALKB8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	493562C3Rik;83431D3Rik;94388N1Rik;ABH8;alkB alkylation repair homolog 8 (E. coli);alkB alkylation repair homolog 8;AlkB homologue 8;ALKB8_HUMAN;ALKBH8;Alkylated DNA repair protein alkB homolog 8;FLJ3824;MGC1235;OTTHUMP211286;OTTHUMP211287;Probable alpha-ketoglutarate-dependent dioxygenase ABH8;RGD134687;S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8;tRNA (carboxymethyluridine(34)-5-O)-methyltransferase ABH8;ALKB8;ALKBH8 ABH8;Alkylated DNA repair protein alkB homolog 8 (EC 1.14.11.-) (Probable alpha-ketoglutarate-dependent dioxygenase ABH8) (S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8) (tRNA (carboxymethyl
Observed Band	
Calculated Molecular Weight	75kD
Cell Pathway	Cytoplasm . Nucleus . Predominantly cytoplasmic. .
Tissue Specificity	Widely expressed, with highest expression in spleen, followed by pancreas and lung.
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images