



ODR4 mouse mAb

Catalog No	YP-mAb-09757
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ODR4 C1orf27 TTG1 TTG1A
Immunogen	Synthesized peptide derived from human ODR4 AA range: 25-75
Specificity	This antibody detects endogenous levels of human ODR4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ODR4;ODR4 C1orf27 TTG1 TTG1A;Protein odr-4 homolog (hODR-4) (LAG1-interacting protein) (Transactivated by transforming growth factor beta protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images