



S41A2 mouse mAb

Catalog No	YP-mAb-09759
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC41A2
Immunogen	Synthesized peptide derived from human S41A2 AA range: 122-172
Specificity	This antibody detects endogenous levels of human S41A2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S41A2;SLC41A2;Solute carrier family 41 member 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images