



UBR4 mouse mAb

Catalog No	YP-mAb-09764
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UBR4 KIAA0462 KIAA1307 RBAF600 ZUBR1
Immunogen	Synthesized peptide derived from human UBR4 AA range: 809-859
Specificity	This antibody detects endogenous levels of human UBR4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	600 kDa retinoblastoma protein-associated factor;E3 ubiquitin-protein ligase UBR4;EC:2.3.2.27;N-recognin-4;p600;UBR4;UBR4 KIAA0462 KIAA1307 RBAF600 ZUBR1;E3 ubiquitin-protein ligase UBR4 (EC 6.3.2.-) (600 kDa retinoblastoma protein-associated factor) (N-recognin-4) (Retinoblastoma-associated factor of 600 kDa) (RBAF600) (p600) (Zinc finger UBR1-type pr
Observed Band	570 kDa, 270 kDa
Calculated Molecular Weight	5183 aa, 570 kDa
Cell Pathway	Membrane ; Multi-pass membrane protein . Cytoplasm. Cytoplasm, cytoskeleton. Nucleus. Concentrates at the leading edge of membrane structures involved in actin motility.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images