



TRGC2 mouse mAb

Catalog No	YP-mAb-09766
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRGC2 TCRGC2
Immunogen	Synthesized peptide derived from human TRGC2 AA range: 127-177
Specificity	This antibody detects endogenous levels of human TRGC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRGC2;TRGC2 TCRGC2;T-cell receptor gamma-2 chain C region (T-cell receptor gamma chain C region PT-gamma-1/2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images