



TF3C1 mouse mAb

Catalog No	YP-mAb-09768
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GTF3C1
Immunogen	Synthesized peptide derived from human TF3C1 AA range: 176-226
Specificity	This antibody detects endogenous levels of human TF3C1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZp686A111;GTF3C1;TF3C alpha;TFIIIC;TFIIIC 220 kDa subunit;TFIIIC box B binding subunit;TFIIIC220;TFIIICalpha;TF3C1;General transcription factor 3C polypeptide 1 (TF3C-alpha) (TFIIIC box B-binding subunit) (Transcription factor IIIC 220 kDa subunit) (TFIIIC 220 kDa subunit) (TFIIIC220) (Transcription factor IIIC su
Observed Band	250 kDa
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images