



TCEA3 mouse mAb

Catalog No	YP-mAb-09769
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TCEA3 TFIISH
Immunogen	Synthesized peptide derived from human TCEA3 AA range: 222-272
Specificity	This antibody detects endogenous levels of human TCEA3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TFIIS;TFIIS.H;TFIISH;Transcription elongation factor S-II protein 3;TCEA3;TCEA3 TFIISH;Transcription elongation factor A protein 3 (Transcription elongation factor S-II protein 3) (Transcription elongation factor TFIIS.h)
Observed Band	39 kDa
Calculated Molecular Weight	348 aa, 39 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images