



MP17L mouse mAb

Catalog No	YP-mAb-09773
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MPV17L
Immunogen	Synthesized peptide derived from human MP17L AA range: 39-89
Specificity	This antibody detects endogenous levels of human MP17L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MP17L;MPV17L;Mpv17-like protein (M-LP homolog) (M-LPH)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 1]: Peroxisome membrane ; Multi-pass membrane protein.
Tissue Specificity	Isoform 1 is detected in the kidney (at protein level). Isoform 1 and isoform 2 are expressed in the kidney, heart, liver, lung, pancreas and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images