



GC224 mouse mAb

Catalog No	YP-mAb-09776
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GCRG224
Immunogen	Synthesized peptide derived from human GC224 AA range: 1-24
Specificity	This antibody detects endogenous levels of human GC224
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AQP 4;AQP-4;AQP4;AQP4_HUMAN;Aquaporin type 4;Aquaporin-4;Aquaporin4;HMIWC 2;HMIWC2;Mercurial insensitive water channel;Mercurial-insensitive water channel;MGC22454;MIWC;WCH 4;WCH4;GC224;GCRG224;Putative gastric cancer-related gene 224 protein (GCRG-224) (GCRG-p224)
Observed Band	
Calculated Molecular Weight	34 kDa; 35kD
Cell Pathway	
Tissue Specificity	Expressed in gastric mucosa.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images