



STOX2 mouse mAb

Catalog No	YP-mAb-09778
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STOX2 KIAA1392
Immunogen	Synthesized peptide derived from human STOX2 AA range: 329-379
Specificity	This antibody detects endogenous levels of human STOX2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1392;Storkhead-box protein 2;STOX2;STOX2 KIAA1392
Observed Band	103 kDa
Calculated Molecular Weight	926 aa, 103 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images