



GARL3 mouse mAb

Catalog No	YP-mAb-09779
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GARNL3
Immunogen	Synthesized peptide derived from human GARL3 AA range: 705-755
Specificity	This antibody detects endogenous levels of human GARL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GARL3;GARNL3;GTPase-activating Rap/Ran-GAP domain-like protein 3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images