



EST1A mouse mAb

Catalog No	YP-mAb-09781
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SMG6 C17orf31 EST1A KIAA0732
Immunogen	Synthesized peptide derived from human EST1A AA range: 688-738
Specificity	This antibody detects endogenous levels of human EST1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C17orf31;EST1-like protein A;Ever shorter telomeres 1A;hSmg5/7a;KIAA732; Smg-6 homolog;SMG6;Telomerase subunit EST1A;Telomerase-binding protein EST1A;EST1A;SMG6 C17orf31 EST1A KIAA0732;Telomerase-binding protein EST1A (EC 3.1.-.-) (EST1-like protein A) (Ever shorter telomeres 1A) (Smg-6 homolog) (Telomerase subunit EST1A) (hSmg5/7a)
Observed Band	
Calculated Molecular Weight	16kDa; 16kD
Cell Pathway	Nucleus, nucleolus . Chromosome, telomere . Cytoplasm, cytosol . Particularly enriched in the nucleolus. .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images