



T200A mouse mAb

Catalog No	YP-mAb-09783
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM200A KIAA1913 HBE61
Immunogen	Synthesized peptide derived from human T200A AA range: 236-286
Specificity	This antibody detects endogenous levels of human T200A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T200A;TMEM200A KIAA1913 HBE61;Transmembrane protein 200A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in cerebellum.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images