



CC137 mouse mAb

Catalog No	YP-mAb-09784
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCDC137
Immunogen	Synthesized peptide derived from human CC137 AA range: 174-224
Specificity	This antibody detects endogenous levels of human CC137
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serumby affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	32kD
Synonyms	CC137;CCDC137;Coiled-coil domain-containing protein 137
Observed Band	
Cell Pathway	Chromosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

