



LUZP6 mouse mAb

Catalog No	YP-mAb-09785
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LUZP6 MPD6
Immunogen	Synthesized peptide derived from human LUZP6 AA range: 8-58
Specificity	This antibody detects endogenous levels of human LUZP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LUZP6;LUZP6 MPD6;Leucine zipper protein 6 (Myeloproliferative disease-associated 6 kDa antigen)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Widely expressed, highest levels found in brain, placenta, spleen, testis, and ovary. Up-regulated in some tumor cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images