



CXB6 mouse mAb

Catalog No	YP-mAb-09787
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GJB6
Immunogen	Synthesized peptide derived from human CXB6 AA range: 118-168
Specificity	This antibody detects endogenous levels of human CXB6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Connexin 3;Connexin-3;Cx3;CXB6_HUMAN;DFNA3;DFNA3B;DFNB1B;ECTD2;ectodermal dysplasia 2;hidrotic (Clouston syndrome);ED2;EDH;Gap junction beta 6 protein;Gap junction beta-6 protein;gap junction protein;beta 6 (connexin 3);beta 6;3kDa;GJB6;HED;HED2;NSRD1;CXB6;Gap junction beta-6 protein (Connexin-30) (Cx30)
Observed Band	
Calculated Molecular Weight	3kDa, 4 kDa; 3kD
Cell Pathway	Cell membrane; Multi-pass membrane protein. Cell junction, gap junction.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

