



FA89B mouse mAb

Catalog No	YP-mAb-09789
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM89B
Immunogen	Synthesized peptide derived from human FA89B AA range: 51-101
Specificity	This antibody detects endogenous levels of human FA89B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA89B;FAM89B;Protein FAM89B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 3]: Cytoplasm . Cell projection, lamellipodium . Co-localizes with CDC42BPA, CDC42BPB and LIMK1 in the lamellipodium. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images