



TM102 mouse mAb

Catalog No	YP-mAb-09795
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM102
Immunogen	Synthesized peptide derived from human TM102 AA range: 275-325
Specificity	This antibody detects endogenous levels of human TM102
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM102;TMEM102;Transmembrane protein 102 (Common beta-chain associated protein) (CBAP)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass membrane protein . Although it is found at the cell surface, the majority of the molecules seems to be located in intracellular compartments.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images