



CH25H mouse mAb

Catalog No	YP-mAb-09798
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CH25H
Immunogen	Synthesized peptide derived from human CH25H AA range: 185-235
Specificity	This antibody detects endogenous levels of human CH25H
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CH25H;Cholesterol 25-hydroxylase (EC 1.14.99.38) (Cholesterol 25-monooxygenase) (h25OH)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images