



NAT6 mouse mAb

Catalog No	YP-mAb-09799
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NAT6 FUS2
Immunogen	Synthesized peptide derived from human NAT6 AA range: 107-157
Specificity	This antibody detects endogenous levels of human NAT6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NAA80;NAT6;N-acetyltransferase 6;HsNAAA80;FUS2;NAT6 FUS2;N-acetyltransferase 6 (EC 2.3.1.-) (Protein fusion-2) (Protein fus-2)
Observed Band	45 kDa
Calculated Molecular Weight	31 kDa
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Strongly expressed in heart and skeletal muscle, followed by brain and pancreas, with weak expression in kidney, liver, and lung and no expression in placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images