



MAD3 mouse mAb

Catalog No	YP-mAb-09800
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MXD3 BHLHC13 MAD3
Immunogen	Synthesized peptide derived from human MAD3 AA range: 96-146
Specificity	This antibody detects endogenous levels of human MAD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IKBA;IκB α;NFKBIA;I kappa B alpha;IkappaBalpha;MAD3;MXD3 BHLHC13 MAD3;Max dimerization protein 3 (Max dimerizer 3) (Class C basic helix-loop-helix protein 13) (bHLHc13) (Max-associated protein 3) (Max-interacting transcriptional repressor MAD3) (Myx)
Observed Band	36 kDa
Calculated Molecular Weight	36 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images