



KHDR2 mouse mAb

Catalog No	YP-mAb-09802
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KHDRBS2 SLM1
Immunogen	Synthesized peptide derived from human KHDR2 AA range: 269-319
Specificity	This antibody detects endogenous levels of human KHDR2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KHDR2;KHDRBS2 SLM1;KH domain-containing;RNA-binding;signal transduction-associated protein 2 (Sam68-like mammalian protein 1) (SLM-1) (hSLM-1); RNA-binding; signal transduction-associated protein 2 (Sam68-like mammalian protein 1) (SLM-1) (hSLM-1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in brain, lung, kidney and small intestine. Weakly expressed in placenta, liver, spleen, thymus, ovary and colon.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images