



STH mouse mAb

Catalog No	YP-mAb-09803
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STH
Immunogen	Synthesized peptide derived from human STH AA range: 52-102
Specificity	This antibody detects endogenous levels of human STH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Allergic rhinitis;ALRH;BHR1;Bronchial hyperresponsiveness 1 (bronchial asthma);IL 13;IL-13;II13;IL13_HUMAN;Interleukin 13;Interleukin-13;MGC116786;MGC116788;MGC116789;NC3;P6;STH;Saitohin
Observed Band	
Calculated Molecular Weight	16kDa; 16kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Highest expression in placenta, muscle, fetal brain, and adult brain, with lower expression in heart, kidney, stomach, testis, and adrenal gland. In the central nervous system, highest expression is in temporal lobe, hypothalamus, medulla and spinal cord, with lower expression in other brain regions.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images