



CECR1 mouse mAb

Catalog No	YP-mAb-09805
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CECR1 ADA2 ADGF IDGFL
Immunogen	Synthesized peptide derived from human CECR1 AA range: 243-293
Specificity	This antibody detects endogenous levels of human CECR1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CECR1;CECR1 ADA2 ADGF IDGFL;Adenosine deaminase CECR1 (EC 3.5.4.4) (Cat eye syndrome critical region protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Detected in blood plasma (at protein level). Widely expressed, with most abundant expression in human adult heart, lung, lymphoblasts, and placenta as well as fetal lung, liver, and kidney. In embryo, expressed in the outflow tract and atrium of the developing heart, the VII/VIII cranial nerve ganglion, and the notochord.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images