



KAZD1 mouse mAb

Catalog No	YP-mAb-09807
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KAZALD1 FKSG28 FKSG40 UNQ2945/PRO21184
Immunogen	Synthesized peptide derived from human KAZD1 AA range: 4-54
Specificity	This antibody detects endogenous levels of human KAZD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KAZD1;KAZALD1 FKSG28 FKSG40 UNQ2945/PRO21184;Kazal-type serine protease inhibitor domain-containing protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images