



ERH mouse mAb

Catalog No	YP-mAb-09810
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ERH
Immunogen	Synthesized peptide derived from human ERH AA range: 10-60
Specificity	This antibody detects endogenous levels of human ERH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DROER;ERH;Enhancer of rudimentary homolog
Observed Band	12 kDa
Calculated Molecular Weight	104 aa, 12 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in all tissues examined.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images