



PITM3 mouse mAb

Catalog No	YP-mAb-09812
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PITPNM3 NIR1
Immunogen	Synthesized peptide derived from human PITM3 AA range: 563-613
Specificity	This antibody detects endogenous levels of human PITM3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PITM3;PITPNM3 NIR1;Membrane-associated phosphatidylinositol transfer protein 3 (Phosphatidylinositol transfer protein;membrane-associated 3) (PITPnm 3) (Pyk2 N-terminal domain-interacting receptor 1) (NIR-1); membrane-associated 3) (PITPnm 3) (Pyk2 N-terminal domain-interacting receptor 1) (NIR-1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endomembrane system ; Peripheral membrane protein .
Tissue Specificity	Detected in brain and spleen, and at low levels in ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images