



WFKN1 mouse mAb

Catalog No	YP-mAb-09813
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WFIKKN1 C16orf12 GASP2 WFIKKN PSEC0040
Immunogen	Synthesized peptide derived from human WFKN1 AA range: 83-133
Specificity	This antibody detects endogenous levels of human WFKN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WFKN1;WFIKKN1 C16orf12 GASP2 WFIKKN PSEC0040;WAP;kazal; immunoglobulin;kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP; follistatin;kunitz and N; kazal; immunoglobulin; kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP; follistatin; kunitz and N
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in pancreas, kidney, liver, placenta, and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images