



RBM4B mouse mAb

Catalog No	YP-mAb-09818
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBM4B RBM30
Immunogen	Synthesized peptide derived from human RBM4B AA range: 83-133
Specificity	This antibody detects endogenous levels of human RBM4B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RBM30;RBM4B;RBM4L;RNA binding motif protein 30;RNA binding motif protein 4B;RNA binding protein 30;RNA binding protein 4B;ZCCHC15;ZCRB3B;RBM4B RBM30;RNA-binding protein 4B (RNA-binding motif protein 30) (RNA-binding motif protein 4B) (RNA-binding protein 30)
Observed Band	36-40 kDa
Calculated Molecular Weight	40 kDa
Cell Pathway	Nucleus. Nucleus, nucleolus.
Tissue Specificity	Expressed in liver and kidney (at protein level). Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images