



TDR10 mouse mAb

Catalog No	YP-mAb-09819
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TDRD10 UNQ9380/PRO34205
Immunogen	Synthesized peptide derived from human TDR10 AA range: 253-303
Specificity	This antibody detects endogenous levels of human TDR10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TDR10;TDRD10 UNQ9380/PRO34205;Tudor domain-containing protein 10
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images