



ARP3C mouse mAb

Catalog No	YP-mAb-09820
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACTR3C ARP11
Immunogen	Synthesized peptide derived from human ARP3C AA range: 86-136
Specificity	This antibody detects endogenous levels of human ARP3C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARP3C;ACTR3C ARP11;Actin-related protein 3C (Actin-related protein 11)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in kidney, stomach, spleen, bone marrow, uterus, testis, placenta, skeletal muscle, mammary gland, lung, fetal liver, and fetal kidney, but not detected in small intestine, brain, and thymus. Expressed in low-metastatic lung adenocarcinoma cells but not in high-metastatic ones.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images