



ME15P mouse mAb

Catalog No	YP-mAb-09821
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	METTL15P1 METT5D2
Immunogen	Synthesized peptide derived from human ME15P AA range: 66-116
Specificity	This antibody detects endogenous levels of human ME15P
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ME15P;METTL15P1 METT5D2;Putative methyltransferase-like protein 15P1 (EC 2.1.1.-) (Methyltransferase 5 domain-containing protein 2) (Methyltransferase-like protein 15 pseudogene 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images