



NANG2 mouse mAb

Catalog No	YP-mAb-09825
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NANOGP1 NANOG2
Immunogen	Synthesized peptide derived from human NANG2 AA range: 136-186
Specificity	This antibody detects endogenous levels of human NANG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NANG2;NANOGP1 NANOG2;Putative homeobox protein NANOG2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images