



PWP2 mouse mAb

Catalog No	YP-mAb-09827
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PWP2 PWP2H
Immunogen	Synthesized peptide derived from human PWP2 AA range: 564-614
Specificity	This antibody detects endogenous levels of human PWP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PWP2H;PWP2;PWP2 PWP2H;Periodic tryptophan protein 2 homolog
Observed Band	102 kDa
Calculated Molecular Weight	919 aa, 102 kDa
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images