



CT118 mouse mAb

Catalog No	YP-mAb-09829
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C20orf118
Immunogen	Synthesized peptide derived from human CT118 AA range: 110-160
Specificity	This antibody detects endogenous levels of human CT118
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	3-keto acyl-CoA synthase ELOVL4;ADMD;Cancer/testis antigen 118;CT118; Elongation of very long chain fatty acids (FEN1/Elo2 SUR4/Elo3 yeast) like 4; Elongation of very long chain fatty acids like 4;Elongation of very long chain fatty acids protein 4;ELOV L4;ELOV4_HUMAN;ELOVL 4;ELOVL4;FLJ17667; FLJ92876;Stargardt disease 3;Stargardt disease 3 autosomal dominant;STGD 2; STGD 3;STGD2;STGD3;C20orf118;Uncharacterized protein C20orf118
Observed Band	
Calculated Molecular Weight	37 KD; 37kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images