



TXLNB mouse mAb

Catalog No	YP-mAb-09830
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TXLNB C6orf198 MDP77
Immunogen	Synthesized peptide derived from human TXLNB AA range: 425-475
Specificity	This antibody detects endogenous levels of human TXLNB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MDP77;LST001;hMDP77;C6orf198;Beta-taxilin;TXLNB;TXLNB C6orf198 MDP77;Beta-taxilin (Muscle-derived protein 77) (hMDP77)
Observed Band	75 kDa
Calculated Molecular Weight	684 aa, 77 kDa
Cell Pathway	
Tissue Specificity	Expressed in skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images