



TSPY4 mouse mAb

Catalog No	YP-mAb-09832
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TSPY4
Immunogen	Synthesized peptide derived from human TSPY4 AA range: 90-140
Specificity	This antibody detects endogenous levels of human TSPY4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CASK interacting nucleosome assembly protein;CDA1;Cell division autoantigen 1;CINAP;CTCL;CTCL tumor antigen se2 4;CTCL-associated antigen se2-4;Cutaneous T cell lymphoma associated antigen se2 4;Cutaneous T cell lymphoma associated tumor antigen se2 4;Cutaneous T-cell lymphoma-associated antigen se2-4;Differentially Expressed Nucleolar TGF beta1 Target;Differentially Expressed Nucleolar TGF beta1 Target protein;Differentially-expressed nucleolar TGF-beta1 target protein;HRIHFB2216;NP79;Nuclear protein of 79 kDa;RP1 29F12.2;SE2 4;SE24;Testis specific protein Y encoded like 2;testis-specific protein Y encoded-like 2;Testis-specific Y-encoded-like protein 2;TSPX;TSPY like 2;TSPY like protein 2;TSPY-like protein 2;TSPYL2;TSYL2_HUMAN;TSPY4;Testis-specific Y-encoded protein 4
Observed Band	
Calculated Molecular Weight	79kDa; 79kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images