



MOGT1 mouse mAb

Catalog No	YP-mAb-09834
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MOGAT1 DC2 DGAT2L1
Immunogen	Synthesized peptide derived from human MOGT1 AA range: 209-259
Specificity	This antibody detects endogenous levels of human MOGT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	613A14Rik;11164N14Rik;2-acylglycerol O-acyltransferase 1;Acyl CoA:monoacylglycerol acyltransferase 1;Acyl-CoA:monoacylglycerol acyltransferase 1;DGAT2L;DGAT2L1;Diacylglycerol acyltransferase 2 like; Diacylglycerol acyltransferase 2-like protein 1;Diacylglycerol O acyltransferase 2 like 1;Diacylglycerol O-acyltransferase candidate 2;hDC2;mDC2;MGAT1; MGC11835;MOGAT1;MOGT1;MOGT1_HUMAN;Monoacylglycerol O acyltransferase 1;Monoacylglycerol O-acyltransferase 1;OTTMUSP25181; MOGAT1 DC2 DGAT2L1;2-acylglycerol O-acyltransferase 1 (EC 2.3.1.22) (Acyl-CoA:monoacylglycerol acyltransferase 1) (MGAT1) (Diacylglycerol O-acyltransferase candidate 2) (hDC2) (Diacylglycerol acyltransferase 2-like prot
Observed Band	
Calculated Molecular Weight	39kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in stomach and liver.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images