



CTGE1 mouse mAb

Catalog No	YP-mAb-09836
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CTAGE1 CTAGE2
Immunogen	Synthesized peptide derived from human CTGE1 AA range: 19-69
Specificity	This antibody detects endogenous levels of human CTGE1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cancer/testis antigen 21.2;CT21.1;CT21.2;CTAGE2;CTAGE;CTAGE 1;CTAGE 2;CTAGE1;Protein cTAGE 2;CTGE1;CTAGE1 CTAGE2;Cutaneous T-cell lymphoma-associated antigen 1 (Protein cTAGE-1) (Cancer/testis antigen 21.1) (CT21.1)
Observed Band	86-120 kDa
Calculated Molecular Weight	754 aa, 89 kDa
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Testis. Expressed in several cutaneous T-cell lymphoma (CTCL) cell lines, in head-neck carcinomas and carcinoma of ovarian tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images