



PLPL3 mouse mAb

Catalog No	YP-mAb-09837
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PNPLA3 ADPN C22orf20
Immunogen	Synthesized peptide derived from human PLPL3 AA range: 123-173
Specificity	This antibody detects endogenous levels of human PLPL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PLPL3;PNPLA3 ADPN C22orf20;Patatin-like phospholipase domain-containing protein 3 (EC 3.1.1.3) (Acylglycerol O-acyltransferase) (EC 2.3.1.-) (Adiponutrin) (Calcium-independent phospholipase A2-epsilon) (iPLA2-epsilon)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type II membrane protein . Lipid droplet .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images