



YAED1 mouse mAb

Catalog No	YP-mAb-09840
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	YAE1D1 C7orf36 GK003
Immunogen	Synthesized peptide derived from human YAED1 AA range: 18-68
Specificity	This antibody detects endogenous levels of human YAED1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	YAED1;YAE1D1 C7orf36 GK003;Yae1 domain-containing protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images